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NH IMMUNOCHROMATO *Listeria* <<Instruction Manual>>

* Read this manual before you start using this test kit

[Intended use]

NH Immunochromato *Listeria* is intended for use in food testing laboratories for qualitative detection of *Listeria* species from various food products. Not intended for diagnostic use.

This kit is intended for use by personnel familiar with the appropriate aseptic techniques for the isolation and identification of *Listeria*.

[Introduction]


Listeria is widely distributed in nature and can be found in soil, rivers, and vegetables. Since it is more resistant to heat, salt, and acid than other bacteria and can grow even in refrigerators¹⁾, it may reach infective dose during distribution and storage of foods.

Outbreaks of listeriosis associated with foods including milk, cheese, vegetables, and meat have been reported around the world. Although no foodborne outbreak has been reported in Japan, it was revealed that the prevalence of *Listeria* in retail foods in Japan is similar to that in Western countries^{2,3)}.

NH Immunochromato *Listeria* is a simple and rapid immunochromatographic test kit for detection of *Listeria* species. The kit does not require skilled personnel or well-equipped laboratory.

[Kit contents]

- | | | |
|----|-------------------------|-------------------------------|
| A: | Test strip | 2 tests × 10 packs (20 tests) |
| B: | Instruction manual | 1 sheet |
| C: | Plastic bag with zipper | 1 piece |

[Principle of assay]**1) Illustration of test strip**

- | | |
|----|-------------------------|
| a. | Sample application zone |
| b. | Reagent zone |
| c. | Nitrocellulose membrane |
| d. | Absorbent pad |
| e. | Test zone |
| f. | Control zone |

Note 1: Do not scratch or touch with bare fingers on the sample application zone (a) or nitrocellulose membrane (c).

2) Principle of assay

When a sample solution is applied to the sample application zone of the test strip, it is wicked through reagent zone. *Listeria* (1) in the sample binds to the colloidal gold-labeled anti-*Listeria* spp. antibody (2) in the reagent zone. The pathogen-antibody complex flows through the nitrocellulose membrane which contains a zone of antibody specific to *Listeria* spp. The immune complex is captured and concentrated in the zone (test zone) (3), eventually forms a reddish purple line.

The membrane also contains a control zone where the gold-labeled anti-*Listeria* spp. antibody is captured. The control line forms in the control zone (4) regardless of the presence of *Listeria* in the sample solution. The control line ensures the test is working correctly.

[Detection limit]

The test using 9 strains of *Listeria* (*L. monocytogenes*, *L. innocua*, *L. ivanovii*, *L. seeligeri*, *L. welshmeri*, and *L. grayi*) revealed that enriched sample including 1×10^5 – 10^7 CFU/mL of *Listeria* will yield a positive result.

Note 2: The lower detection limit of the kit may vary depending on *Listeria* serotypes, strains, and components of test solution.

[Precautions]

- (1) Contamination with enriched cultures can cause *Listeria* infection because they can contain high amounts of *Listeria*. Users of NH Immunochromato *Listeria* should be familiar with the appropriate aseptic techniques for the isolation and identification of *Listeria*. Extreme care should be taken in handling samples, enriched cultures and used test strips.
- (2) Do not use expired test strips. The expiration date is indicated on the aluminum pouches and the box.
- (3) Good microbiological laboratory practices should be used to prevent contamination of samples, lab workers and the environment.


[Additionally required materials and instruments]

1. Enrichment broth (e.g. half Fraser broth, Fraser broth, EB broth, UVM broth)
2. Stomacher and stomacher bags (preferably with a filter)
3. Autoclave
4. Incubator (30°C, 37°C)
5. Balance (capable of weighing 25 g)
6. Disposable plastic transfer pipettes and/or appropriate micro pipettes and disposable tips
7. Timer
8. Thermostatic bath (that for which temperature can be set at 80°C or higher)
9. Glass test tubes

[Test procedures]

- (1) Weigh 25-g sample into the stomacher bag and add 225-mL the primary enrichment broth (half Fraser broth).
- (2) Homogenize with a stomacher. The length of time for homogenization may vary depending on sample type.
- (3) Incubate the homogenized sample at 30°C for 24±3 hours.
- (4) Remove the stomacher bag from the incubator and mix the culture gently using side-to-side motion.
- (5) Inoculate 0.1 mL of the primary enrichment culture to 10 mL of the secondary enrichment broth (Fraser broth) and incubate at 37°C for 48±3 hours.
- (6) Remove the test tube from the incubator and mix the culture gently.
- (7) Dispense 5 mL of the enriched culture to a glass test tube with a sterilized pipette. Incubate the test tube at 80°C for 20 minutes. The incubated culture should be applied to test strips after cooling to room temperature.
- (8) Allow the test strips warm to room temperature (20–25°C) before opening the aluminum pouch. Open the pouch immediately before use (Notes 4, 5).

*Enriched sample should be applied to a test strip by either of following procedures (9a or 9b).


- (9a) Place a test strip on flat surface and with a new pipette tip apply 100 µL of the enriched sample, drop-wise to the sample application zone. Careful not to overload the sample application zone; overloading it may cause spillage.
- (9b) Transfer 150 µL of the enriched sample into a test tube with a new pipette tip. Put a test strip in the test tube immersing the sample application zone in the sample solution (Note 6).
- (10) Observe the test result 15 minutes after applying the sample solution to the test strip.

Note 3: Keep the remaining enriched samples unsterilized until the end of the test. Confirmatory tests, such as isolation of *Listeria*, might be required after testing with the kit. Unsterilized enriched samples can be used in confirmatory tests.

Note 4: Do not open the aluminum pouch before it comes to the room temperature to prevent moisture absorption by test strips. Moist test strips can give incorrect test results. Unused test strips should be kept at 2–8°C in the plastic bag (provided with the kit) along with a desiccant sheet.

Note 5: Do not scratch or touch with bare fingers on the sample application zone or nitrocellulose membrane. Hold the absorbent pad when handling test strips.

Note 6: Keep the surface of sample solution below the arrows printed on the test strip.

[Interpretation of the test results]

- (1) Positive results: Reddish purple lines appear in the test zone and the control zone in 15 minutes after application of sample solutions.
- (2) Negative results: A reddish purple line appears in the control zone.
- (3) Invalid results: No reddish purple line appears in the control zone.

Note 7: Presumptive positive results should be confirmed according to authorized methods⁴⁾. The enriched samples prepared for testing with the kit can be used in confirmatory tests.

[Disposal]

Decontaminate used test strips, test tubes, pipettes, and enriched culture by autoclave, etc. in strict compliance with local, state and federal disposal regulations.

[Storage and expiration date]

- (1) Storage: NH immunochromato *Listeria* should be stored at 2–8°C under protection from the light. NH immunochromato *Listeria* must not be frozen.
- (2) Expiration date: 12 months from the date of manufacture.

[Warranty]

NH Foods Ltd. makes no warranty, whether expressed or implied. If any materials are defective, Distributor or NH Foods Ltd. will provide a replacement of the product. There is no warranty of merchantability or fitness for a particular purpose. Each user is expected to verify the appropriateness of each operative procedure and the application to individual food materials. NH Foods Ltd. shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

[References]

- (1) Mitsuyama M.: "Listeria", . Microbiology, Asakura Publishing, 454-460 (2002)
- (2) Ministry of Health, Labor, and Welfare eds.: Standard Methods of Analysis in Food Safety Regulation Food Sanitation Inspection Guidelines (Microbiology), Japanese Association of Food Sanitation Hygiene Association, 249-265 (2004)
- (3) Igimi S.: Food Sanitation Research, 53, 19-23 (2003)
- (4) MHLW Notice No.1128, Article 3, from the Department of Food Safety "Detection Method for *Listeria monocytogenes*"

[Technical assistance]

For technical assistance, please contact distributor or manufacturer.

● **Distributor** (Labeled on the box)

● **Manufacturer**

R&D Center, NH Foods Ltd.

3-3 Midorigahara, Tsukuba, Ibaraki 300-2646, Japan

Phone: +81-29-847-7825 Fax: +81-29-848-1256

E-mail: nh.kit@nipponham.co.jp

http://www.rdc.nipponham.co.jp